


# Pronestander

# LECKEY<sup>®</sup>



## How to unpack and assemble

When you receive your Leckey Pronestander it will be boxed and packed flat.

If there are two people present, it is possible to lift the standing frame out of the top of the box.

If there is only one person available to remove the Pronestander from its box, we would recommend that you:

Open the top of the box to determine which way the frame is sitting.

Carefully use a sharp knife to cut down the two back corner edges of the box.

Unlock the two back swivel castors and wheel the Pronestander out of the box.

### ACCESSORIES:


Chest Plate  
Knee Plate  
Pelvic Plate and Belt  
Hinged Foot Plate and Sandals  
Activity Tray, Poles and Bowl  
Castors  
Angle Clamp Plate/s

## How to assemble

There is a choice of five different ANGLE POSITIONS on the Pronestander. Depending on the user's tolerance to standing, you can choose to position them anywhere from a very prone angle to almost upright.

In order to fix the Pronestander at the desired angle, you should;


1. Lock all the castors, facing out from the frame, at an angle of 45 degrees.



2. Hold the chest plate and gently swing the main frame of the stander upright.



3. Loosen the three safety bar knobs (or six with the adult size which has two clamp plates) sufficiently for you to be able to remove the safety clamp plate/s.



4. Position the angle adjuster bar at the desired angle between the location bars.



5. Slip the safety clamp plate/s back into position and tighten the safety bar knobs securely.




**Never attempt to adjust the angle of the Pronestander when the frame is in use.**

Once you have selected the angle desired, you can set up the rest of the Leckey Pronestander to suit your child. You should first measure the distance from the user's feet to just under their arms. This is best done in a lying position on a mat.



## How to adjust

You can now proceed to set up the standing frame so that the distance between the sandals and the top of the chest plate is equal to this height. This will always be the first adjustment you need to carry out.




To adjust the height of the CHEST PLATE, loosen the 4 main frame knobs (2 on each side of the frame), and this will allow the whole INNER SECTION (including chest plate, pelvic plate and knee plate) to move up and down. Simply select the height required and re-tighten the knobs securely.

N.B. There are only 2 knobs on a size 1 Pronestander.

**Again, never attempt to adjust the height of the main frame when the stander is in use.**

You can then go on to adjust the PELVIC PLATE and KNEE PLATE to the appropriate height for your child's measurements. The pelvic plate is adjusted by loosening the 2 pelvic plate knobs (1 on each side of the plate) and sliding the padded plate into the required position. The knee plate is adjusted in exactly the same way using the 2 knee plate knobs.



The KNEE CUPS can be moved from left to right and each knee cup is fixed by a single knob at the back of the knee plate. Again, loosen this bolt to adjust the knee cup, select the exact position you desire, and re-tighten the knobs securely in place.




Minor adjustments to the knee and pelvic plate may be made during use.

The Pronestander FOOT PLATE is split and hinged at the back of the foot board. This allows the sandals to be adjusted more easily and can facilitate transfer from a seating position, a wheelchair for example.




Each SANDAL is attached to the individual foot plate with velcro and a single fixing knob. If you wish to rotate the sandal laterally, simply loosen this knob and place it in the position you require. The velcro will hold the sandal securely until you re-fasten the fixing knob.

The Pronestander LATERAL SUPPORTS offer three way adjustment. By loosening the 4 lateral support knobs (2 on each side), you can move the supports in and out, and also rotate them in place to produce an angle.



Using the lateral support angle knobs (1 on each side) we can achieve an even finer degree of positioning. You have the option of applying a little more pressure from behind, by loosening these knobs, bringing the supports around the back of your child's trunk, and fixing them securely in place.



The PELVIC BELT should be attached to your child first, again preferably on a mat. The pelvic straps are attached firmly around your child's pelvis, and then, when the user is secured in the frame (by their feet and chest), the longer frame straps are slipped through the D-RINGS on the pelvic plate and secured firmly behind the user with the velcro on the belt.




To attach the TRAY, click the end of the tray poles into the white tray pole brackets positioned at either side of the location bars. There are also 2 of these tray pole brackets on the base of the tray and they attach in exactly the same way.

The tray can then be attached to the body of the Pronestander by slotting the tray pins into any one of the 3 tray holes on each side of the main frame.

Simply pull the 2 tray pins out, away from the tray, get the tray into position and release the pins which will spring back into position and secure the tray.

To increase or decrease the height of the front of the tray (and thereby alter the tray surface angle), loosen the height adjustment knobs at the front of the tray poles and select the height you require.



Your frame comes complete with a TRAY COVER and BOWL, meaning that you can use both the tray with a flat surface, or unclick the cover and use the bowl for play activities.



## How to maintain

We would recommend that users are never transferred straight from seating to standing without preparation and stretching. This is best done in lying whilst on a mat, but can be carried out in a wheelchair or seating system if required.

Lock the castors with the brake lever facing out from the base. Please do not use your hands.

We would recommend that the castor brakes be locked at all times when the frame is in use, being adjusted or even just stored.

Although mobility is an important feature of the Pronestander, we do not recommend that children be moved over uneven surfaces whilst in the equipment.

Never adjust either the angle of the Pronestander, or the overall height of the chest plate whilst the frame is in use. Other finer positional adjustments may be carried out safely when the child is in the frame.

When placing your child into the standing frame, for safety reasons, always secure your child's feet and the chest straps first.

We recommend that children should not be left unattended at any time whilst using Leckey equipment.

At least once a week, check that all fixings and fastenings are in good condition.

The Pronestander is manufactured from lacquered plywood and upholstered in wipe clean vinyl. When cleaning the frame we would recommend that you use warm water and a non-abrasive detergent.

When Pronestander parts are not in use, we would advise that you store them **together** in a safe place for future use.

If you are in any doubt as to the safe use of your Leckey Pronestander, or you require further assistance, please do not hesitate to contact our Customer Service department.



### JAMES LECKEY DESIGN INC

360 MERRIMACK STREET  
RIVER WALK, BUILDING 9  
LAWRENCE, MA 01843

TEL : (508) 691-5220  
FAX : (508) 691-5230

Tollfree : USA 1 800 532-5390  
Freephone : UK 0800 318285  
Web Site : [www.leckey.com](http://www.leckey.com)

### JAMES LECKEY DESIGN LTD

KILWEE BUSINESS PARK,  
DUNMURRY BT17 0HD,  
NORTHERN IRELAND

TEL : 01232 602277  
FAX : 01232 600795

Ireland 1800 626060  
Email : [info@leckey.com](mailto:info@leckey.com)